


α-AMYLASE

(CNP-G3) Liquid Reagent



Intended Use : Kit for the quantitative determination of α-Amylase in human serum.

Contents: Pack Size : 1 x 10 ml

Reagent-1 10 Test

Amylase Reagent : 1 x 10 ml

Insert : 01 No.

Directions for Use:
Use the process instruction outlined in the kit insert

Precautions:

1. For *in-vitro* diagnostics use only .
2. Do not inhale the reagent or pipette by mouth.
3. Avoid contact with skin/eyes. In case of contact, wash off immediately with plenty of water and seek medical attention as early as possible.
4. Do not use if package is damaged.
5. Dispose waste product in accordance with local Govt. regulations.

Storage and Stability: 2°C - 8°C

The unopened reagent is stable up to the expiry date mentioned on the label when stored between 2 - 8°C.


MLC.No. : S(0079)/15/MB

LOT :

M.R.P. Rs. : (incl. of all taxes)


Manufactured by:
Vanguard Diagnostics Private Limited
C-123, Phase-I, Okhla Industrial Area
New Delhi-110 020, INDIA. Tel:011-4172 7222
Email:info@vanguarddiagnostics.com

PM3000BLS2



α-AMYLASE

(CNP-G3)
Liquid Reagent
(Reagent-1)



Ready to Use

STORE AT 2 - 8°C. 2°C - 8°C

Manufactured by:
Vanguard Diagnostics Private Limited
C-123, Phase-I, Okhla Industrial Area
New Delhi-110 020, INDIA.
Tel:011-4172 7222
Email:info@vanguarddiagnostics.com

For *in vitro* Diagnostic use only .

MLC.No. : S(0079)/15/MB

LOT :

M.R.P. Rs. : (incl. of all taxes)

Manufactured by:
Vanguard Diagnostics Private Limited
C-123, Phase-I, Okhla Industrial Area
New Delhi-110 020, INDIA. Tel:011-4172 7222
Email:info@vanguarddiagnostics.com

PM3000BLS3
