



Dengue Combo (NS1/IgG/IgM) Antigen & Antibody Detection Rapid Card Test


Pack Size: 1 x 10 Tests

Intended Use: One step immunochromatographic test for detection of NS1 Antigen and IgG & IgM Antibodies of Dengue Virus in Human Serum/Plasma.

Directions for Use:

Follow the process instructions outlined in the Product Insert.

Precautions:

1. For *in-vitro* diagnostics use only. 
2. Do not use if package is damaged.
3. Dispose waste product in accordance with local Govt. regulations.

Contents:

Test Device	: 10 Nos.
Assay Buffer for Dengue IgG/IgM	: 01 No.
Dropper-1 for Dengue NS1	: 10 Nos.
Dropper-2 for Dengue IgG/IgM	: 10 Nos.
Product Insert	: 01 No.



Vanguard Diagnostics Private Limited
C-123, Phase-I, Okhla Industrial Area
New Delhi-110 020, INDIA.
Tel: 011-4172 7222, Email: info@vanguarddiagnostics.com



PHD0000195