



## HDL/LDL Calibrator



For *in vitro* Diagnostic use only .

MLC.No. : MFG/IVD/2021/000002

STORE AT 2° - 8°C.



Manufactured by:  
Vanguard Diagnostics Private Limited  
C-123, Phase-I, Okhla Industrial Area  
New Delhi-110 020, INDIA.

PIM30008L977

Reconstitute with 1.0 ml Distilled water

HDL Conc.-

LDL Conc.-



:



:



:

M.R.P. ₹ :

(incl. of all taxes)