



LDL CHOLESTEROL

(Direct Reagent Kit) Liquid Reagent



Intended Use : Kit for the quantitative determination of LDL Cholesterol in human serum and plasma.

Contents:

- Reagent 1 : R1 Reagent : 2 x 30 ml
- Reagent 2: R2 Reagent : 2 x 10 ml
- Insert : 01 No.

Pack Size : 2 x 40 ml
▽ : 160 Test


Storage and Stability: 2°C - 8°C

The unopened reagent is stable up to the expiry date mentioned on the label when stored between 2 - 8°C.

Directions for Use:

Use the process instruction outlined in the kit insert

Precautions:

1. For *in-vitro* diagnostics use only .
2. Do not inhale the reagent or pipette by mouth.
3. Avoid contact with skin/eyes. In case of contact, wash off immediately with plenty of water and seek medical attention as early as possible.
4. Do not use if package is damaged.
5. Dispose waste product in accordance with local Govt. regulations.

MLC.No. : S(0079)/15/MB

 :

 :

 :

M.R.P. Rs.:
(Incl. of all taxes)

Manufactured by:
Vanguard Diagnostics Private Limited
C-123, Phase-I, Okhla Industrial Area
New Delhi-110 020, INDIA. Tel:011-4172 7222
Email:info@vanguarddiagnostics.com

PM30008L63



LDL CHOLESTEROL

(Direct Reagent Kit) Liquid Reagent
(Reagent 1)



Ready to use

For *in vitro* Diagnostic use only .

STORE AT 2 - 8°C. 2°C - 8°C

MLC.No. : S(0079)/15/MB

Manufactured by:
Vanguard Diagnostics Private Limited
C-123, Phase-I, Okhla Industrial Area
New Delhi-110 020, INDIA. Tel:011-4172 7222
Email:info@vanguarddiagnostics.com

PM30008L64
 :
 :
 :

30 ml



LDL CHOLESTEROL

(Direct Reagent Kit) Liquid Reagent
(Reagent 2)



Ready to Use

For *in vitro* Diagnostic use only .

STORE AT 2 - 8°C. 2°C - 8°C

MLC.No. : S(0079)/15/MB

Manufactured by:
Vanguard Diagnostics Private Limited
C-123, Phase-I, Okhla Industrial Area
New Delhi-110 020, INDIA.
Tel:011-4172 7222
Email:info@vanguarddiagnostics.com

PM30008L65
 :
 :
 :

10 ml