



Lyse - G

Pack Size: 1 x 500 L

Intended Use: Lyse - G is used for lysing of cells for measurement of WBC & Hemoglobin in Genius, 3-part differential hematology analyzers.

Active Ingredients:

- Sodium Chloride: 1.2%
- Potassium Chloride: 0.03%
- Quaternary Ammonium salt: 3.0%
- Preservative: 0.01%

Directions for Use:

Open the container cover. Pull the connector of the container out and connect it with the corresponding Lyse tubing as per instructions given in the Analyzer's user manual.

Precautions:

1. For *in-vitro* diagnostics use only.
2. Do not inhale the reagent or pipette by mouth.
3. Avoid contact with skin/eyes. In case of contact, wash off immediately with plenty of water and seek medical attention as early as possible.
4. Do not use if package is damaged.
5. Dispose waste product in accordance with local Govt. regulations.

Storage and Stability:

The unopened reagent is stable up to the expiry date mentioned on the label when stored between 15°C - 30°C. Use Lyse - G, within 6 months of first opening.

MLC.No. : MFG/IVD/2021/000002



:



:



:

M.R.P (Rs.) :

(Incl. of all Taxes)

Manufactured by

 Vanguard Diagnostics Private Limited

C-123, Phase-I, Okhla Industrial Area, New Delhi-110020

Tel: 011-4172 7222, Email: info@vanguarddiagnostics.com



PM30008L654