

**VDx™**

# Probe Cleaner - M

**Intended Use:** Probe Cleaner – M is used for periodic cleaning of the tubes & baths of Mindray 3-part differential hematology analyzers.

**Pack Size:** 1x100 ML

**Active Ingredients:**

- o Sodium Hydroxide: 0.5%
- o Sodium hypochlorite solution: 2.5%
- o Hydrochloric Acid: 0.2%

**Directions for Use:**

Open the container cover of Probe Cleaner – M during the system shutdown, proceed as per instructions given in the user manual of the Analyzer.

**Precautions:**

1. For *in-vitro* diagnostics use only.
2. Do not inhale the reagent or pipette by mouth.
3. Avoid contact with skin/eyes. In case of contact, wash off immediately with plenty of water and seek medical attention as early as possible.
4. Do not use if package is damaged.
5. Dispose waste product in accordance with local Govt. regulations.
6. Cap the bottle tightly after use.

**Storage and Stability:**

The unopened reagent is stable up to the expiry date mentioned on the label when stored between 15°C - 30°C. Use Probe Cleaner - M, within 6 months of first opening.

**MLC.No. : MFG/VD/2021/00002**



:



:



:

**M.R.P. (Rs.) :**



:

**Wingard Diagnostics Private Limited**  
100, Shree Chhatrapati Shivaji Maharaj  
Road, Vashi, Mumbai - 401302  
Tel: 022-41533333, Email: [info@wingarddiagnostics.com](mailto:info@wingarddiagnostics.com)



PM3000BL372