



URIC ACID

(URICASE-POD) Liquid Reagent



Intended Use : Kit for the quantitative determination of Uric Acid in human serum and plasma.

Contents:

Reagent-1

Uric Acid Reagent : 4 x 50 ml

Reagent-2

Standard : 1 x 3 ml

Insert : 01 No.

Directions for Use:

Use the process instruction outlined in the kit insert

Precautions:

1. For *in-vitro* diagnostics use only **IVD**.
2. Do not inhale the reagent or pipette by mouth.
3. Avoid contact with skin/eyes. In case of contact, wash off immediately with plenty of water and seek medical attention as early as possible.
4. Do not use if package is damaged.
5. Dispose waste product in accordance with local Govt. regulations.

Pack Size : 4 x 50 ml

▽ : 200 Test

Storage and Stability:



The unopened reagent is stable up to the expiry date mentioned on the label when stored between 2 - 8°C.

MLC.No. : S(0079)/15/MB

LOT :

:

:

M.R.P. :
(Incl. of all taxes)

Manufactured in India by:
Vanguard Diagnostics Private Limited
C-123, Phase-I, Okhla Industrial Area
New Delhi-110 020, INDIA. Tel:011-4172 7222
Email:info@vanguarddiagnostics.com

PM30008L33



URIC ACID

(URICASE-POD) Liquid Reagent

(Reagent 1)



Ready to use

For *in vitro* Diagnostic use only **IVD**.

50 ml

STORE AT 2 - 8°C.

MLC.No. : S(0079)/15/MB

Manufactured in India by:
Vanguard Diagnostics Private Limited
C-123, Phase-I, Okhla Industrial Area
New Delhi-110 020, INDIA. Tel:011-4172 7222
Email:info@vanguarddiagnostics.com

PM30008L34 **LOT** :
 :
 :



URIC ACID

(URICASE-POD) Liquid Reagent

(Reagent 2)



STANDARD

6 mg/dl

Ready to use

STORE AT 2 - 8°C.

For *in vitro* Diagnostic use only **IVD**.

MLC.No. : S(0079)/15/MB

3 ml

Manufactured in India by:
Vanguard Diagnostics Private Limited
C-123, Phase-I, Okhla Industrial Area
New Delhi-110 020, INDIA. Tel:011-4172 7222
Email:info@vanguarddiagnostics.com

PM30008L35 **LOT** :
 :
 :